

QuReader® 8 Green

<http://de.lumiprobe.com/p/fluorometer-qureader-8-well-green>

The QuReader® 8 Green fluorometer is a desk fluorometer with a touch screen designed for dsDNA, ssDNA, RNA, and protein quantitation. The QuReader® fluorometer enables the fluorescence measuring for up to eight samples using green channel. Total compatibility with the QuDye® assays ensure highly sensitive determination of nucleic acid and protein concentrations in the lower pg/uL range.

Allgemeine Eigenschaften

Rechtliche Hinweise: Dieses Produkt wird nur für Forschungszwecke angeboten und verkauft. Es wurde nicht auf Sicherheit und Wirksamkeit in Nahrungsmitteln, pharmazeutischen Produkten, medizinischen Vorrichtungen, Kosmetika sowie für gewerbliche oder andere Einsatzzwecke getestet. Der Verkauf gewährt oder impliziert nicht die Erlaubnis zur Verwendung in der In-vitro-Diagnostik, bei der Herstellung von Nahrungsmitteln oder pharmazeutischen Produkten, in medizinischen Vorrichtungen sowie in kosmetischen Erzeugnissen.