

Magnetic Beads (100 mg/mL)

<http://de.lumiprobe.com/p/magnetic-beads>

Magnetic silica-based particles represent the modern standard for rapid and efficient DNA and RNA isolation. This technology delivers preparations of the highest purity, free from inhibitors and with minimal sample loss. Offering seamless compatibility with automated platforms, it is equally suited for manual protocols and high-throughput research. You will obtain pure, ready-to-use nucleic acids, ideal for even the most demanding downstream applications.

Allgemeine Eigenschaften

Rechtliche Hinweise: Dieses Produkt wird nur für Forschungszwecke angeboten und verkauft. Es wurde nicht auf Sicherheit und Wirksamkeit in Nahrungsmitteln, pharmazeutischen Produkten, medizinischen Vorrichtungen, Kosmetika sowie für gewerbliche oder andere Einsatzzwecke getestet. Der Verkauf gewährt oder impliziert nicht die Erlaubnis zur Verwendung in der In-vitro-Diagnostik, bei der Herstellung von Nahrungsmitteln oder pharmazeutischen Produkten, in medizinischen Vorrichtungen sowie in kosmetischen Erzeugnissen.